

RASTRUM

Protocol | Cell recovery from hydrogel matrices

Document Number | WET.BCM.0009.PRO

Version Number | v1.02

Date of issue | 23rd June 2020

INVENTIA
INSPIRING SCIENCE

Introduction

RASTRUM Cell Recovery Solution is an enzymatic-based reagent for use in rapid extraction of 3D cell models encapsulated in RASTRUM hydrogel for use in downstream applications.

Materials

Kit storage and Handling

- Store reagents at -20 °C upon receipt.
- For best results, do not leave reagents at room temperature (RT) for longer than 1 hr.

Materials required but not provided

- RASTRUM 3D cell models
- 1 x phosphate-buffered saline (PBS) solution

Protocol

1. Thaw reagents, 20-30 min at RT.
2. Aspirate and discard media from printed 3D cell models.
3. Wash with 150 µL 1 x PBS, 5 min at 37 °C.
4. Aspirate and discard PBS.
5. Overlay 75 µL Cell Recovery Solution (Support Fluid F175).
6. Incubate 10 min at 37 °C.

Note: At this stage, dissolution of the biofunctional matrix will be complete, but the inert base will remain intact.




Figure 1. MCF-7 spheroids encapsulated in 3D bioprinted hydrogels. (A) Encapsulated spheroids before dissociation, with the dashed line indicating the hydrogel matrix. (B) Spheroids remain intact after a 10 minute incubation with RASTRUM Cell Retrieval Reagent. The dotted line indicates the location of the former hydrogel matrix. (C) Empty wells post-cell retrieval show efficient hydrogel dissociation and cell retrieval. (D) Previously encapsulated spheroids remained intact after dissociation from the hydrogel matrix. Scale bars = 250 µm.

7. Transfer solution containing dissociated cells into a 1.5 mL eppendorf tube.

Optional: To maximise cell recovery, wash each well with an additional 150 µL PBS and combine with cell suspension from Step 7

8. Pellet cells using cell-specific centrifugation parameters.
9. Remove and discard supernatant.
10. Resuspend cell pellet in appropriate reagent for downstream analysis.

INVENTIA

Inventia Life Science Operations Pty Ltd
ABN 19 613 078 710
Suite 1.13, 90-96 Bourke Road,
Alexandria, NSW, 2015, Australia
Telephone +61 (0) 412 175 725
Fax +61 (2) 8399 03 04
info@inventia.life | www.inventia.life

The products are intended only for laboratory research purposes. They are not to be used for any other purposes, including but not limited to in vitro diagnostic purposes, in foods, drugs, medical devices or cosmetics for humans or animals, or for commercial purposes. The customer warrants that it will not use the products for any such purpose. For further inquiries, please contact technical service.